

SUTURA	IFSA DIVISION ATRAMAT	J&J DIVISION ETHICON	SYNETURE	B. BRAUN
ACIDO POLIGLICOLICO	PGA	VICRYL	DEXON POLYSORB	SAFIL
ACIDO POLIGLICÓLICO DE ABSORCIÓN RÁPIDA	PGA RAPID	VICRYL RAPID	-----	SAFIL QUICK
POLIDIOXANONA	PDX	PDS II	MAXON	MONOPLUS
CATGUT CROMICO	CATGUT CROMICO	CATGUT CROMICO	CATGUT CROMICO	-----
CATGUT SIMPLE	CATGUT SIMPLE	CATGUT SIMPLE	CATGUT SIMPLE	-----
NYLON	NYLON	ETHILON	DERMALON	DAFILON
POLIESTER (SIN RECUBRIMIENTO)	POLIESTER	MERSILENE	SURGIDAC	DAGROFIL
POLIESTER (CON RECUBRIMIENTO)	POLIESTER	ETHIBOND EXCEL	TI-CRON	SYNTHOFIL
POLIPROPILENO	POLIPROPILENO	PROLENE	SURGI PRO	PREMILENE
SEDA (PRECORTADA)	SEDA	SUTUPACK	SOFSILK	SILKAM
SEDA (CON AGUJA)	SEDA	PERMAHAND	SOFSILK	SILKAM
ACERO	ACERO	ACIFLEX	STEEL	STEELEX
POLIGLECAPRONE -25	PGC-25	MONOCRYL	-----	MONOSYN
FLUORURO DE POLIVINILIDENO	PVDF	PRONOVA	-----	-----

EQUIVALENCIA DE AGUJAS

Punta Ahusada 1/2 círculo			
Atramat	J+J	D+G	
F10	RB-3	T-37	
F13	RB-2	T-30	
F17	RB1	T-31	
R20	SH-1	T-16	
G22	CT-3	T-3	
R26	SH	T-5	
R30	MH-1	V-30	
R37	MH	T-10	
G37	CT-1	T-12	
G40	CT	T-25	
G48	CTX	T-20	

Punta Ahusada 3/8 círculo			
Atramat	J+J	D+G	
FE12	C-1	TE-1	
RE30			

Punta Ahusada 5/8 círculo			
Atramat	J+J	D+G	
GH37	UR-5	TT-3	

Recta Ahusada			
Atramat	J+J	D+G	
T50	ST	ST-1	

Reverso Cortante Giratoria 1/2 círculo			
Atramat	J+J	D+G	
J48			
J57			

Reverso Cortante 1/2 círculo			
Atramat	J+J	D+G	
C12	M-1	PR-1	
C13		PR-2	
C15	PS-4		
C18	PS-4	PR-4	
C20	J-1	C-4	
C26	CP-2	C-6	
C30	CP-2	C-30	
C37	CP-1	C-9	
C40	CP	C-23	
C48	CPX	C-12	
C57	CCS-1		

Punta Roma 1/2 círculo			
Atramat	J+J	D+G	
Z65	BP-1	BSG-27	

Reverso Cortante 1/2 círculo			
Atramat	J+J	D+G	
CE12	C-3	CE-2	
CE16	FS-3	CE-3	
CE19	FS-2	CE-4	
CE24	FS-1	CE-6	
CE26	FS	C-15	
CE30	FSL	CE-10	
CE37	FSLX	CE-12	
CE45			
CE77	LR	CE-24	

Cuerpo Redondo Punta Cortante 1/2 círculo - circle			
Atramat	J+J	D+G	
W40			
W48	V40	C-12	
W57			

Recta Cortante			
Atramat	J+J	D+G	
U60	KS	CS-1	

Electrodo para Marcapasos Temporal			
Atramat	J+J	D+G	
UL60	SKS	CS-1	
UL90	SKS3	CS-6	

Cortante Convencional 3/8 círculo			
Atramat	J+J	D+G	
NE19	PC-5	SCE-4	
NE24			

Cuerpo Redondo Punta Cortante Recta			
Atramat	J+J	D+G	
TW254			

Ahusada Loop 1/2 círculo			
Atramat	J+J	D+G	
X40			
X50	CTX		
X65	TP-1		

AGUJAS PREMIUM

Punta Ahusada Premium 1/2 círculo			
Atramat	J+J	D+G	
S10	BV-1	CV-337	
S13	RB-2	CV-330	
S16	RB1	CV-331	
S20	SH-1	CV-316	
S25	CT-2	T-6	
S30	RH plus		

Ahusada Punta Cortante Premium 1/2 círculo			
Atramat	J+J	D+G	
Y16	V-5	KV-16	
Y18			
Y25	V-7	Y-5	

Ahusada Punta Cortante Premium 3/8 círculo			
Atramat	J+J	D+G	
YE13	CC-1	KV-11	
YE16			
YE37			

Cortante Convencional Premium 3/8 círculo			
Atramat	J+J	D+G	
QE13	PC-1		
QE16	PC-3		
QE19	PC-5		
QE20	PC-5	SCE-4	
QE24			
QE26			

Reverso Cortante Premium 1/2 círculo			
Atramat	J+J	D+G	
P7	P-2	RR-1	
P16	PS-4	P24	

Reverso Cortante Premium 3/8 círculo			
Atramat	J+J	D+G	
PE13	P-3	PRE-2	
PE16	PS-3	PRE-3	
PE19	PS-2	PRE-4	
PE24	PS-1		
PE26	PS		

Punta Ahusada Atralloy 3/8 círculo			
Atramat	J+J	D+G	
SBE10	BV	CV-1	
SBE13	C-1	CV-11	

Ahusada Regular Premium 1/2 círculo			
Atramat	J+J	D+G	
SR13			
SR17			
SR20			
SR26			

Ahusada Gruesa Premium 1/2 círculo			
Atramat	J+J	D+G	
SG30			
SG36			
SG40	CT plus		

Punta Ahusada Atralloy (Antirreflejante Negra) 3/8 círculo			
Atramat	J+J	D+G	
SVE8			
SVE10	CN	CV-1	
SVE13	CN-1		

Espatulada Premium 1/4 círculo			
Atramat	J+J	D+G	
EA6	TG-8	SE-100-8	
EA7	TG-100-8		
EA8	RD-1	DO-1	
EA9	S-14	SS-14	

Espatulada Premium 1/2 círculo			
Atramat	J+J	D+G	
E6	TG-160-6	SE-160-8	
E7	TG-175-8	SE-175-8	

Espatulada Premium 3/8 círculo			
Atramat	J+J	D+G	
EE6	TG140-6	LE-140-6	
EE6	TG140-6	LE-140-8	
EE7	TG-140-8	LE-140-8	
EE8			

Punta Ahusada Premium 3/8 círculo			
Atramat	J+J	D+G	
SE4	BV75-3	TE-70	
SE5	BV130-4	TE-100	
SE6	BV130-5	TE-145	
SE7			
SE8	BV175-6	MV-175-8	
SE10	BV	TE-11	
SE13	C-1	CV-11	